

Technical Data Sheet

Product

**Anti-Ms CD8a Purified**

Cat. Number/Size

**11-579-C025**

**0.025 mg**

**11-579-C100**

**0.1 mg**

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

---

Antigen	CD8a
Clone	53-6.7
Format	Purified
Reactivity	Mouse
Application	FC (QC tested), IP, IHC(P), IHC(F), ICC
Application details	Flow cytometry: Recommended dilution: 1.5 µg/ml. Immunohistochemistry (frozen sections): Recommended dilution: 1:1000; formaldehyde fixation is not recommended, acetone fixation is preferred.
Isotype	Rat IgG2a kappa
Specificity	The rat monoclonal antibody 53-6.7 recognizes an extracellular epitope of murine CD8a (32-34 kDa; alpha chain of the CD8 antigen).
Other names	Ly-2, Ly-B, Ly-35, Lyt-2, BB154331
Immunogen	Murine spleen cells
Entrez Gene ID	12525
Gene name	CD8a
NCBI Full Gene Name	CD8a molecule
UniProt ID	P01731

---

Concentration	1 mg/ml
Preparation	Purified by protein-G affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.